

HeartSine® samaritan® PAD 350P/360P AEDs

Philips HeartStart OnSite



HeartSine samaritan PAD 350P/360P

Philips HeartStart OnSite

✓	Compact, portable, lightweight 1.1 kg (2.4 lbs)	✗	Larger, heavier 1.5 kg (3.3 lbs)
✓	Highest level of protection from dust and water ingress in the industry (IP56)	✗	Lower water and dust resistance rating (IP21)
✓	Integrated battery and electrodes	✗	Separate battery and electrodes
✓	4-year electrode life/4-year battery life	✗	2-year electrode life/4-year battery life
✓	One item (Pad-Pak™) to track and replace after use or upon expiration	✗	Two items to track and replace after use or upon expiration
✓	Fully automatic version (SAM 360P)	✗	No fully automatic version
✓	90-minute ECG recording time	✗	15-minute ECG recording time
✓	Low energy escalating protocol (SCOPE™)	✗	Fixed non-escalating energy protocol
✓	Operating impedance range of 20 – 230 ohms	✗	Operating impedance range of 25 –180 ohms
✓	Audible and visual prompts to help guide rescuer in a noisy environment	✗	Audible prompts only
✓	Carry case included	✗	Carry case not included
✓	Field upgradeable software	✗	Not field upgradeable
✓	Free Saver EVO™ data management software (data cable required)	✗	No free data management software (data cable required)
✓	Paediatric compatible (When used with Paediatric-Pak™)	✓	Paediatric compatible
✓	Automatic self test	✓	Automatic self test
✓	8-year warranty	✓	8-year unit warranty
✓	Lowest total cost of ownership; acquisition and operating expense		

EMEA/APAC

HeartSine Technologies, Ltd.

203 Airport Road West
Belfast, Northern Ireland
BT3 9ED

Tel: +44 28 9093 9400

Fax: +44 28 9093 9401

info@heartsine.com



U.S./Americas

HeartSine Technologies LLC

121 Friends Lane, Suite 400
Newtown, PA 18940

Toll Free: (866) 478 7463

Tel: +1 215 860 8100

Fax: +1 215 860 8192

info@heartsine.com



The HeartSine products described in this brochure meet the European Medical Directive requirement.



UL Classified.
See complete marking on product.

H009-032-247-1

All names herein are trademarks or registered trademarks of their respective owners.
© 2017 HeartSine Technologies LLC. All rights reserved.

www.heartsine.com



HeartSine®